

2017 CDT Codes Installation Guide for CS PracticeWorks

Installing CDT 2017 Codes to the CS PracticeWorks Software

This job aid contains instructions for installing the 2017 CDT code set to the CS PracticeWorks software, version 7.10.0 or higher.

CS PracticeWorks v7.10.0 contains the ADA's new CDT 2017 code set.

Note: To request the CS PracticeWorks v7.10.0 update, go to www.carestreamdental.com/updatecentral.

To install the 11 new 2017 codes into the 2017 CDT code set in v7.10.0:

- 1 On a computer running CS PracticeWorks, v7.10.0 or higher, select **Start > All Programs > CS PracticeWorks > Utilities**. The **CS PracticeWorks Utilities** window is displayed.
- 2 Select **PATCHES**. The **Patch Master** window is displayed.
- 3 Select the patch labeled **CDT 2017, Install - Update Transaction Code List with CDT 2017 codes**, and click **Run selected patch**. An End User License Agreement Window is displayed.
- 4 **Select I understand and agree to this license agreement**, and click **Agree**. The **Install CDT 2017 Codes** window is displayed with a list.
- 5 Click **Yes**. A progress screen is displayed, followed by an **Installation Complete** message.
- 6 Click **OK**.
- 7 Return to the Patch Master window and select the patch labeled **CDT 2017, Enable - Enable CDT 2017 codes for all insurance companies**, and click **Run selected patch**. This patch enables all insurance companies to accept the new code set without displaying "Are you sure?" messages.
- 8 Close the **Patch Master** and **Utilities** windows.

To confirm that the update was successful:

- 1 Open the PracticeWorks Software and select **Lists > Transaction Codes**. The **Transaction Code List** is displayed.
- 2 Click the **Code #** sorting option at the top of the second column.
- 3 In the **Find** field, type **0414**. This code denotes **Laboratory Processing of Microbial Specimen to Include Culture and Sensitivity Studies, Preparation and Transmission of Written Report**, which is the first of the new CDT 2017 codes.

If you have issues during the installation, contact a Carestream Dental Support representative at 800.603.4438, option 1.

